EXHIBIT 9

FiberWire™

IMPORTANT PRODUCT INFORMATION WICHTIGE PRODUKTINFORMATION NOTICE D'UTILISATION IMPORTANTE IMPORTANTI INFORMAZIONI PER L'USO INSTRUCCIONES IMPORTANTES PARA EL USO

CE 0086

Manufacturer:
Arthrex Inc.
Naples, Florida 34108-1945 • USA
Toll-Free: +1 800 934-4404
www.arthrex.com

E. C. Representative:
Arthrex Med. Inst. GmbH
85757 Karlstfeld
Germany

Tel: +49 81 31 59 57 0 • Fax: +49 81 31 59 57 63 1

DFU-0085
Rev. 5

Description:

Arthrex FiberWire is available in several U.S.P. sizes (sutures meet U.S.P. standards for suture, except diameter). The Arthrex FiberWire may also be sold with needles attached (swaged) to the ends in a variety of sizes. The suture is made of polyethylene fibers and polyester fibers braided, sterilized and coated for surgical use. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue. The Arthrex FiberWire is available non-dyed (white) or dyed and meets or exceeds U.S.P. and European standards (except for diameter).

Indications:

Arthrex FiberWire is indicated for use in soft tissue approximation and/or ligation. FiberWire is not for use in cardiac indications.

Actions:

Arthrex FiberWire, when tested per ISO/DIS 10993, Biological Evaluation of Medical Devices-Part 10: Tests for Irritations and Sensitization, had no reactions of allergic or sensitive nature. The dyed suture and coating are pharmacologically inactive.

Arthrex FiberWire is not absorbed, but may become encapsulated in the surrounding connective tissues. The Arthrex FiberWire is not known to have significant change in tensile strength in vivo.

Contraindications:

None known

Warnings:

Do not re-sterilize. Once open, discard unused suture. Do not expose to heat.

Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Arthrex FiberWire for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

As with any foreign body, prolonged contact of this or any other suture with soft solutions, such as those found in the urinary or biliary tracts, may result in calculus formations. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Precautions:

In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Assure that all knots have been secured using accepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Care should be taken to prevent damage to surrounding

tissue or user puncture due to improper handling of the needlepoint.

Do not grasp the needle at the point or swage, to avoid damage to these areas. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in "sharps" containers.

Adverse Reactions:

Adverse reactions have not been noted with the Arthrex FiberWire product in animal testing. Common non-absorbable suture reactions may include wound dehiscence, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, enhanced bacterial infectivity, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

Sterilization:

Arthrex FiberWire suture is supplied sterile. Method of sterilization - EO. Do not resterilize. Do not use if package is opened or damaged. Discard opened, unused sutures.


Storage Conditions:

Store below 25°C, away from moisture and direct heat. Do not use after expiration date.

How Supplied:

The Arthrex FiberWire is available in several U.S.P. sizes (sutures meet U.S.P. standards for suture, except diameter). The suture is supplied sterile in pre-cut lengths and in some cases with swaged needles. The Arthrex FiberWire is available in non-dyed (white) or dyed colors. The suture is made of polyethylene fibers and polyester fibers braided, sterilized and coated for surgical use. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue.

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