

BROOKSTEIN DECLARATION EXHIBIT 15

CE
0085



FiberWire™

IMPORTANT PRODUCT INFORMATION WICHTIGE PRODUKTINFORMATION NOTICE D'UTILISATION IMPORTANTE IMPORTANTI INFORMAZIONI PER L'USO INSTRUCCIONES IMPORTANTES PARA EL USO

Manufacturer:
Arthrex Inc.
Naples, Florida 34108-1945 • USA
Tel/Fax: +1 800 334-4404
www.arthrex.com

E. C. Representative:
Arthrex Med. Inst. GmbH
65757 Karlshof
Germany
Tel: +49 81 31 59 57 0 • Fax: +49 81 31 59 57 63 1



DU-0085
Rev. 3

Description:
Arthrex FiberWire is available in several U.S.P. sizes (sutures meet U.S.P. standards for suture, except diameter). The Arthrex FiberWire may also be sold with needles attached (designed to be used in a variety of sizes). The suture is made of polyethylene fibers and polyester fiber braided, sutured and coated for surgical use. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue. The Arthrex FiberWire is available non-dyed (white) or dyed in shades of colors (U.S.P. and European standards (except for diameter)).

Indications:
Arthrex FiberWire is indicated for use in soft tissue approximation and/or ligation. FiberWire is not for use in cardiac indications.

Actions:
Arthrex FiberWire, when tested per ISO 10993, Biological Evaluation of Medical Devices - Part 12: Tests for Irritation and Sensitization, had no reactions of allergic or sensitive nature. The dyed suture and coating are pharmacologically inactive.

Arthrex FiberWire is not absorbed, but may become encapsulated in the surrounding connective tissues. The Arthrex FiberWire is not known to have significant change in tensile strength in vivo.

Contraindications:
None known

Warnings:
Do not re-sterilize. Once open, discard unused suture. Do not expose to heat.

Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing Arthrex FiberWire for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used.

As with any foreign body, prolonged contact of this or any other suture with soft tissues, such as those found in the urinary or biliary tracts, may result in calcific formations. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Precautions:
In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

Assure that all ends have been secured using accepted surgical knot tying techniques. Adequate knot security requires the accepted surgical technique of full, square ties, with additional throws as warranted by surgical circumstances and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments. Care should be taken to prevent damage to surrounding

Issue or user procedure due to improper handling of the needlepoint.

Do not grasp the needle at the point or attempt to exert damage to these areas. Handling needles may cause them to lose strength and be less resistant to bending and breaking. Discard used needles in "sharp" containers.

Adverse Reactions:

Adverse reactions have not been noted with the Arthrex FiberWire product in animal testing. Common non-absorbable suture reactions may include wound dehiscence, calcific formation in urinary and biliary tracts when prolonged contact with soft solutions such as urine and bile occurs, enhanced bacterial infection, minimal acute inflammatory tissue reaction, pain, edema, and erythema at the wound site. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of bloodborne pathogens.

Sterilization:

Arthrex FiberWire suture is supplied sterile.

Method of sterilization: EO

Do not resterilize. Do not use if package is opened or damaged. Discard opened, unused sutures.


Storage Conditions:

Store below 25°C, away from moisture and direct heat. Do not use after expiration date.

How Supplied:

The Arthrex FiberWire is available in several U.S.P. sizes (sutures meet U.S.P. standards for suture except diameter). The suture is supplied sterile in pre-cut lengths and in some cases with tapered needles. The Arthrex FiberWire is available in non-dyed (white) or dyed colors. The suture is made of polyethylene fibers and polyester fibers braided, sterilized and coated for surgical use. The coating acts as a lubricant for suture sliding, knot tying, and ease of passing suture through tissue.

SYMBOLS USED ON LABELING



Beschreibung:

Arthrex FiberWire ist in verschiedenen USP-Größen erhältlich (das Normenmaterial entspricht dem USP-Normen für Nahtmaterial, mit Ausnahme des Durchmessers). Arthrex FiberWire ist unter Umständen auch mit an der Fadenenden versehenen (geschliffenen) Nadeln in unterschiedlichen Größen erhältlich. Das Nahtmaterial besteht aus gewickelten, sterilisierten und für den chirurgischen Gebrauch bestimmten Polyethylen- und Polyesterfasern. Die Beschichtung fungiert als Gleitmittel für das Durchdringen von Gewebe und erleichtert das Durchziehen der Fäden durch das Gewebe. Arthrex FiberWire ist ungerichtet (weiß) oder gefärbt. Arthrex FiberWire ist sterilisiert (EO) und entspricht den europäischen Standards (mit Ausnahme des Durchmessers).

Anwendungsgebiete:

Arthrex FiberWire ist für Weichteilapproximation und/oder -ligation vorgesehen. FiberWire nicht für kardio-indikationen verwenden.

Funktionen:

Tests bei Arthrex FiberWire gemäß ISO 10993, Biologischer Evaluation of Medical Devices - Part 12: Prüf- und Sensibilisierungstests ergaben keine allergischen oder empfindlichen Reaktionen. Das gefärbte Nahtmaterial und die Beschichtung sind pharmakologisch inaktiv.

Arthrex FiberWire wird zwar nicht absorbiert, jedoch kann Umkapselung von umgebenden Bindegewebsgeweben. Bei Arthrex FiberWire wurde in vivo keine signifikante Änderung der Zugsfestigkeit festgestellt.

Gegenanzeigen:

Unbekannt

Warnhinweise:
Nach sterilen, unbenutzten Fadenmaterial nicht den Öfen entsorgen. Vor Hitze warnen.

Benutzer sollten vor dem Verschieben von Wunden mit Arthrex FiberWire mit den chirurgischen Präparations-Techniken vertraut sein. Bei denen nicht-überwachten Fäden verwendet wird, ist das Dehiscenzrisiko je nach Anwendungsstelle und verwendetem Fadenmaterial unterschiedlich zu sein.

Wie bei Fremdörpern aller Art kann der längere Kontakt dieses oder jedes anderen Fadenmaterials mit Schilddrüse, wie sie z.B. im Harn- und Gallenkanal vorhanden sind, zu Calcifizierung führen. Bei der Drainage und beim Schließen von Infektionen oder Kontaminationen Wunden sind die in der Chirurgie üblichen Praktiken zu beachten.

Vorsichtsmaßnahmen:

Bei der Handhabung dieses oder jedes anderen Fadenmaterials sorgfältig darauf achten, dass das Material nicht beschädigt wird. Schützen durch Zusammenklappen oder Abklemmen mit chirurgischen Instrumenten wie Zangen oder Nadeln nach Möglichkeit vermeiden.

Sicherstellen, dass sämtliche Knoten gemäß den akzeptierten chirurgischen Knotenbindungspraktiken sicher befestigt wurden. Voraussetzung für angemessene Knotensicherheit ist die Verwendung von Zeugen, qualitativ hochwertigen Nahtmaterialien und Einwirkung des Zugspannung.

Chirurgien. Besonders beim Verfügen von monofilen Fäden ist eine Umkapselung des Fadenendes zu vermeiden. Bei Arthrex FiberWire wurde in vivo keine signifikante Änderung der Zugsfestigkeit festgestellt.

Das Nahtmaterial ist in verschiedenen USP-Größen erhältlich. Das Nahtmaterial besteht aus gewickelten, sterilisierten und für den chirurgischen Gebrauch bestimmten Polyethylen- und Polyesterfasern. Die Beschichtung fungiert als Gleitmittel für das Durchdringen von Gewebe und erleichtert das Durchziehen der Fäden durch das Gewebe. Arthrex FiberWire ist ungerichtet (weiß) oder gefärbt. Arthrex FiberWire ist sterilisiert (EO) und entspricht den europäischen Standards (mit Ausnahme des Durchmessers).

Nebenwirkungen:

Bei Tierversuchen wurden bei der Verwendung von Arthrex FiberWire keine Nebenwirkungen festgestellt. Zu den bei nicht-absorbierbaren Fäden üblichen Reaktionen zählen unter Umständen Drüsenentzündungen, Calcifizierung in Harn- und Gallenwegen bei längerem Kontakt mit Sekretionsorganen (wie sie im Harn und in der Gallenblase vorhanden sind), verstärkte Bakterieninfektion, minimale akute Gewebsentzündungen, Schmerzen, Ödem und Erythem an der Wundstelle. Vereinzelt können Stiche mit kontaminierten chirurgischen Nadeln kann zur Übertragung von Blutpathogenen führen.

Stabilität:

Arthrex FiberWire wird nicht absorbiert, jedoch kann Umkapselung von umgebenden Bindegewebsgeweben. Bei Arthrex FiberWire wurde in vivo keine signifikante Änderung der Zugsfestigkeit festgestellt.


Lagerungsbedingungen:

Unter 25°C trocknen und fern von direkter Hitze einwirkung lagern. Nicht nach dem Verfallsdatum verwenden.

Literatur:

Arthrex FiberWire ist in verschiedenen USP-Größen erhältlich. Das Nahtmaterial besteht aus gewickelten, sterilisierten und für den chirurgischen Gebrauch bestimmten Polyethylen- und Polyesterfasern. Die Beschichtung fungiert als Gleitmittel für das Durchdringen von Gewebe und erleichtert das Durchziehen der Fäden durch das Gewebe.

AUF DER VERPACKUNG VERWENDETE SYMBOLE



BROOKSTEIN DECLARATION

EXHIBIT 16

TO : Arthrex

ATTN: Don Grafton

FROM : Brian Hallett

DATE : 19/10/2000

SUBJECT : Polyester - Dyneema Braid

Dear Don,

Please find enclosed 4 DT trials samples for your inspection ,these have been made using Polyester/Dyneema mixed either in the cover or straight core, to match US2 I have set out below a matrix of how each was made and their results for your information

DT PA23 SAMPLE COMMENTS: CORE DID NOT BREAK ON KNOT PULL ONLY COVER, STRAIGHT PULL CORE BROKE COVER STAYED INTACT.

16 Carrier m/c

COVER 16 carriers in use each with 1 end of 138 d'tex Polyester per carrier

CORE 1 end of 165/1/3 with 10 TPI "S" and 7 TPI "Z"(Dyneema)

PPI 38

Stage	St/pull kg	Knot/pull kg	Runnage mt/kg	Diameter mm	Extensin %	Solids %
M/c	19.12	9.87	3455	0.677	5.48	
Dye	17.05	8.81		0.589	10	
Stretch	16.51	6.95		0.577	5.3	15.3
Finish	17.2	10.35	3375	0.569	6.8	

DEPUY MITEK
EXHIBIT 164
04cv12457

PR 06515
Confidential - Outside Attorneys' Eyes Only

DT PA25 SAMPLE COMMENTS: BOTH CORE AND COVER BROKE ON THE STRAIGHT PULL AT DYE STAGE, ON THE STRAIGHT PULL AT STRETCH STAGE THE CORE BROKE ONLY

16 Carrier m/c

COVER 16 carriers in use

8 car with 1 end of 113 poly 8 car with 1 end of 110 dyneema

CORE 1 end of 190/1/3 with 10 TPI "S" and 7 TPI "Z"

PPI 50

Stage	St/pull kg	Knot/pull kg	Runnage m/kg	Diameter mm	Extensin%	Solids %
M/c	28.82	11.293	3803	0.582	11.15	
Dye	25.56	10.5		0.582	14.7	
Stretch at 5%	26.84	10.58		0.587	9.9	15.3
Finish	24.35	11.95	3703	0.55	11.2	

DT PA26 SAMPLE COMMENTS: CORE BROKE IN FIRST TWO READINGS, WHOLE BRAID BROKE IN THIRD

16 Carrier m/c

COVER 16 carriers in use

8 car with 1 end of 113 polyester 8 car with 1 end of 110 Dyneema

CORE 1 end of 165/1/3 with 10 TPI "S" and 7 TPI "Z"(Dyneema)

PPI 50

Stage	St/pull kg	Knot/pull kg	Runnage m/kg	Diameter mm	Extensin%	Solids %
M/c	21.82	10.953	3908	0.681	5.28	
Dye	23.62	12.26		0.693	10.1	
Stretch at 5%	24.56	11.79		0.573	5.8	15.3
Finish	21.48	12.87	3786	0.578	6	

DT PA27 SAMPLE COMMENTS: CORE DID NOT BREAK ON KNOT PULL ONLY COVER, STRAIGHT PULL CORE BROKE COVER STAYED INTACT.

16 Carrier m/c

COVER 16 carriers in use

16 car with 1 end of 113 Polyester per carrier

CORE 1 end of 165/1/3 with 10 TPI "S" and 7 TPI "Z"(Dyneema)

PPI 44

Stage	St/pull kg	Knot/pull kg	Runnage m/kg	Diameter mm	Extensin%	Solids %
M/c	19.14	9.037	4033	0.549	5.9	
Dye	16.18	8.35		0.548	8.5	
Stretch at 5%	16.49	6.54		0.545	5.5	15.3
Finish	16.12	8.04	3919	0.553	5.6	

If I can be any further assistance please do not hesitate to contact me

Kind regards

Brian Hallett

Brian Hallett

Product Development Manager

PR 06514

Confidential - Outside Attorneys' Eyes Only

BROOKSTEIN DECLARATION EXHIBIT 17

IN THE UNITED STATES DISTRICT COURT
FOR THE DISTRICT OF MASSACHUSETTS

DePuy Mitek, Inc., a
Massachusetts Corporation,

Plaintiff,

vs.

CIVIL ACTION
NO. 04-12457 PBS

Arthrex, Inc., a Delaware
Corporation,

Defendant.

DEPOSITION OF:

ASHLEY HOLLOWAY

DATE:

September 15, 2005

TIME:

1:08 p.m. to 5:07 p.m.

LOCATION:

The Ritz Carlton Golf Resort
2600 Tiburon Drive
Naples, FL 34112

TAKEN BY:

Plaintiff

REPORTER:

Deborah A. Krotz, RPR, CRR

VIDEOGRAPHER:

Les Smoak, CLVS

<p>1 suture --</p> <p>2 A. Correct.</p> <p>3 Q. -- with respect to knot security; is that right?</p> <p>4 A. Correct.</p> <p>5 Q. Okay. Does Arthrex test its FiberWire sutures</p> <p>6 for pliability?</p> <p>7 A. No.</p> <p>8 Q. Do you know what pliability means?</p> <p>9 A. Pliability --</p> <p>10 Q. As it relates to FiberWire sutures?</p> <p>11 A. I don't know the exact definition, no.</p> <p>12 Q. But Arthrex does not test its FiberWire sutures</p> <p>13 for pliability?</p> <p>14 A. No.</p> <p>15 Q. Has it ever tested its FiberWire sutures for</p> <p>16 pliability?</p> <p>17 A. Not that I'm aware of.</p> <p>18 Q. Does Arthrex test its FiberWire sutures for</p> <p>19 handleability?</p> <p>20 A. Yes.</p> <p>21 Q. How does Arthrex test its FiberWire sutures for</p> <p>22 handleability?</p> <p>23 A. It's a subjective test.</p> <p>24 Q. What do you mean by that?</p> <p>25 A. I mean basically we give a piece of suture to a</p>	<p>30</p> <p>1 Q. What do you mean by that?</p> <p>2 A. I mean everything that's in that construct.</p> <p>3 Q. Contributes to the handleability of the suture?</p> <p>4 A. Yes.</p> <p>5 Q. What is the handleability of Arthrex's FiberWire</p> <p>6 suture?</p> <p>7 MR. TAMBURIO: Objection to form.</p> <p>8 A. Do you want me to give you a subjective answer?</p> <p>9 Q. Let me rephrase the question. Has Arthrex</p> <p>10 received feedback from surgeons on the handleability of</p> <p>11 Arthrex's FiberWire suture?</p> <p>12 A. Yes.</p> <p>13 Q. And what was the feedback that Arthrex received</p> <p>14 from surgeons on the handleability of Arthrex's FiberWire</p> <p>15 suture?</p> <p>16 A. That it's easy to utilize.</p> <p>17 Q. And what surgeons provided that feedback?</p> <p>18 A. I have heard Dr. Burkhardt say that.</p> <p>19 Q. Anyone else?</p> <p>20 A. Not directly to me, no.</p> <p>21 Q. Have any doctors provided negative feedback to</p> <p>22 Arthrex on the handleability of Arthrex's FiberWire</p> <p>23 suture?</p> <p>24 A. Yes.</p> <p>25 Q. How many doctors have provided negative feedback</p>
<p>31</p> <p>1 product manager or a surgeon who is familiar with the</p> <p>2 field of sutures and ask them to give us feedback on the</p> <p>3 handleability.</p> <p>4 Q. And what is handleability measured in?</p> <p>5 A. There are no units. It's subjective.</p> <p>6 Q. So what does "handleability" mean? What does the</p> <p>7 definition of "handleability" mean as Arthrex uses that as</p> <p>8 a test for its FiberWire sutures?</p> <p>9 A. What do I think it means?</p> <p>10 Q. What does Arthrex think it means?</p> <p>11 A. I think handleability as is it easy to move</p> <p>12 through the tissue or pass through the tissue? Is it easy</p> <p>13 to slide knots? Is it easy to tie knots?</p> <p>14 Q. Anything else?</p> <p>15 A. Is it easy to slide through the anchor eyelet?</p> <p>16 Q. Anything else?</p> <p>17 A. Not that I can think of.</p> <p>18 Q. Do materials contribute to the handleability of</p> <p>19 Arthrex's FiberWire sutures?</p> <p>20 MR. TAMBURIO: Objection to form.</p> <p>21 A. Yes.</p> <p>22 MR. TAMBURIO: Also, it's outside the scope.</p> <p>23 Q. What materials contribute to the handleability of</p> <p>24 Arthrex's FiberWire sutures?</p> <p>25 A. All materials used.</p>	<p>33</p> <p>1 to Arthrex on the handleability of Arthrex's FiberWire</p> <p>2 suture?</p> <p>3 A. I know of only one account.</p> <p>4 Q. And who was that?</p> <p>5 A. I believe, again, it might have been Dr.</p> <p>6 Burkhardt.</p> <p>7 Q. And what negative feedback did Dr. Burkhardt</p> <p>8 provide to Arthrex on the handleability of Arthrex's</p> <p>9 FiberWire sutures?</p> <p>10 A. That the nylon was too repetitious in the</p> <p>11 TigerWire.</p> <p>12 Q. What does that mean?</p> <p>13 A. It means there were too many wraps.</p> <p>14 Q. And that affected the handleability of the</p> <p>15 suture?</p> <p>16 A. That affected the feel.</p> <p>17 Q. Is the feel like handleability, or is that</p> <p>18 another word for handleability -- feel?</p> <p>19 A. Feel would fall under handleability; how does it</p> <p>20 feel in your hands.</p> <p>21 Q. All right. So that's one criteria used for the</p> <p>22 handleability?</p> <p>23 A. Right.</p> <p>24 Q. As well as passing the suture through tissue or</p> <p>25 any other things you described?</p>

BROOKSTEIN DECLARATION EXHIBIT 18

1 IN THE UNITED STATES DISTRICT COURT
2 FOR THE MIDDLE DISTRICT OF MASSACHUSETTS

3 DEPUY MITEK, INC.,
4 a Massachusetts corporation,

5 Plaintiff,

6 v.

Case No: CA-0412457-PBS

7 ARTHREX, INC, a
8 Delaware corporation,

9 Defendant.
10 _____ /

11 VIDEOTAPE DEPOSITION OF ANN WATERHOUSE

12 TAKEN: Pursuant to Notice by
13 Counsel for the Plaintiff

14 PLACE: Ritz Carlton Golf Resort
15 2600 Tiburon Drive
16 Naples, FL 34109

17 DATE: Wednesday, August 24, 2005

18 TIME: Began: 8:55 a.m.
19 Ended: 1:00 p.m.

20 BEFORE: TRACIE L. MOUNTAIN-THOMPSON
21 Court Reporter
22 Notary Public
23 State of Florida at Large
24
25

COPY

Page 46

1 Q What does that mean?

2 A In the first paragraph, under suture
3 weight, which is coating, it says, "Pearsalls makes the
4 following statement about the NuSil coating applied to
5 both our Arthrex FiberWIRE and other competitive
6 product."

7 So they're demonstrating their processing of
8 that coating by the statement that follows.

9 Q That statement is from NuSil, not from Arthrex?

10 A No. That's from Pearsalls.

11 Q All right. But just to be clear, the coating
12 on the FiberWIRE product is MED-2174, right?

13 A Correct.

14 Q Okay. And in that last paragraph on page
15 2104, the second sentence says, "As noted (above) we
16 cannot measure the amount of coating so the product is
17 accepted by our customers on the basis of an agreed
18 detailed coating process which includes mixing the NuSil
19 to a certain viscosity and the speeds, temperatures, and
20 other parameters for the coating process. For each
21 coating batch all these details are recorded in the batch
22 documentation available to you and other customers."

23 Do you see that?

24 A Yes.

25 Q What does that mean?

Page 48

1 A Correct.

2 Q "We," Pearsalls, "cannot measure the amount of
3 coating so the product is accepted by our customers on
4 the basis of an agreed detailed coating process which
5 includes mixing the NuSil to a certain viscosity and the
6 speeds, temperatures and other parameters for the coating
7 process."

8 A Correct.

9 Q So does Arthrex accept batches of FiberWIRE
10 based on the coating?

11 MR. SABER: Objection. Vague. Confusing
12 question.

13 BY MR. FALKE:

14 Q All right. Was Pearsalls representing to
15 Arthrex that it cannot measure the amount of coating on
16 the FiberWIRE product when Arthrex submitted this
17 statement to the FDA in Exhibit 81?

18 A Can you repeat the question.

19 Q Sure.

20 Was Pearsalls representing to Arthrex that
21 Pearsalls cannot measure the amount of coating on the
22 FiberWIRE product when Arthrex submitted Exhibit 81 to
23 the FDA?

24 A Yes, as a percentage.

25 Q As a percentage of what?

Page 47

1 A That means that that's the processing that
2 they -- they put onto our product and it's describing
3 that they do it to a certain viscosity. They do it at a
4 certain speed and temperature. And they have other
5 parameters for the coating process. So that coating
6 process is an agreed upon process by their buyers,
7 basically, so that's what it is. That's how it's
8 actually done.

9 Q Does Arthrex and Pearsalls have an agreed -- an
10 agreed detailed coating process?

11 A Not that I know of.

12 Q So how does Arthrex know how the FiberWIRE
13 coating is applied to its FiberWIRE product?

14 A Specifically, they could ask for that batch
15 record for each lot that is produced, but in
16 general, they know the coating process because of the
17 temperature and speed, and the parameters that Pearsalls
18 uses coats that in an even manner, so it's an assumption.

19 Q Is it true that Pearsalls cannot measure the
20 amount of coating on the FiberWIRE product?

21 MR. SABER: Objection. Inconsistent with the
22 testimony.

23 BY MR. FALKE:

24 Q I'm just looking. It says -- on page 2104 it
25 says, "We," and I assume that's Pearsalls; is that right?

Page 49

1 A The total weight of the suture.

2 Q Has Arthrex submitted any document to the FDA
3 in any submission which details the coating process
4 Pearsalls uses to coat any FiberWIRE suture?

5 A No.

6 Q Why not?

7 A Because we used the statement as a reference
8 point so that we didn't have to submit anything about the
9 coating process. We described it instead.

10 Q You say "the statement," which statement is
11 that?

12 A On 2104, the statement that Pearsalls makes.

13 Q Okay. And then it says, in the last sentence
14 on page 2104, "For each coating batch all these details
15 are recorded in the batch documentation available to you
16 and other customers."

17 Do you see that?

18 A Yes.

19 Q Is the reference to "you" in that sentence,
20 does that refer to Arthrex?

21 A Correct.

22 Q Does Arthrex have any documents which reflect
23 the details of the coating process for the batches of
24 FiberWIRE it receives?

25 A Not that I know of.

13 (Pages 46 to 49)

1 Q Does Arthrex have the ability to obtain those
2 documents?

3 A Yes.

4 Q So Pearsall's can determine the percentage --
5 cannot determine the percentage of coating to the total
6 weight of the suture; is that right?

7 A Correct.

8 MR. SABER: Objection. Vague.

9 BY MR. FALKE:

10 Q I'm sorry, what was your answer?

11 A Correct.

12 Q And the reason they can't do that is because
13 they can't measure the amount of coating on it;
14 therefore, they can't measure the percentage of the
15 amount of coating to the total weight, right?

16 A Correct.

17 Q I'm going to reask that question because the
18 record is slightly unclear. So Pearsall's cannot
19 determine the percentage of coating to the total weight
20 of the suture; is that right?

21 MR. SABER: Could you read that question back,
22 please.

23 (The court reporter read back the requested
24 portion of testimony.)

25 MR. SABER: Objection. Vague. Asked and

1 answered.

2 BY MR. FALKE:

3 Q You can answer.

4 A Correct.

5 Q Okay. If you could turn to page 2130, please,
6 in Exhibit 81.

7 I'm sorry, are you there?

8 A Uh-huh.

9 Q If you look at the last paragraph on that page,
10 2130, and in particular, the last sentence, it says --
11 excuse me, let me step back. The second to the last
12 paragraph on page 2130, it says, in the last sentence of
13 that paragraph, "The dyed polyester suture (D & C Blue
14 No. 6) and di-peroxide silicone oil (coating) are
15 pharmacologically inactive."

16 Do you see that?

17 A Yes, I do.

18 Q Does the di-peroxide silicone oil (coating)
19 refer to the MED-2174 that's applied to the FiberWIRE
20 product?

21 A Yes, it does.

22 Q And what does that mean "pharmacologically
23 inactive"?

24 A Pharmacologically inactive means that it does
25 not cause any reaction when tested.

1 Q What do you mean by "when tested"?

2 A We test according to 10993, which is cited in
3 that paragraph.

4 Q Uh-huh.

5 A There are certain tests that you have to
6 perform according to that based on the use of your
7 product. When they test that, as an example -- and
8 you've got both of the tests -- they test whether there
9 is a reaction or -- a negative reaction or an adverse
10 reaction as opposed to nothing happening.

11 Q Right.

12 A If it -- if there is nothing happening, they
13 consider that pharmacologically inactive.

14 Q Has Arthrex made any representation to the FDA
15 in any submission as to whether MED-2174 is absorbable in
16 the body?

17 A No.

18 Q Do you know if MED-2174 is absorbable in the
19 body?

20 A No.

21 Q No, you don't know?

22 A I don't know.

23 Q I think I asked this but I'm going to ask it
24 again. And Exhibit 81 was submitted to the FDA by
25 Arthrex, right?

1 A Yes, it was.

2 (Exhibit 82 was marked for identification.)

3 BY MR. FALKE:

4 Q Let me hand you Exhibit 82, which is a document
5 with Bates numbers ARM001882 through 1884. Have you seen
6 this document before?

7 A Yes, I have.

8 Q And what is Exhibit 82?

9 A This is the substantial equivalence letter for
10 K010673.

11 Q And that reference, that K number, 010673,
12 reverts to Exhibit 78, right?

13 A Correct.

14 Q And what is the purpose of this substantial
15 equivalence letter, Exhibit 82?

16 A It's the FDA's way of granting you permission
17 for sale of a product.

18 Q Okay. So through Exhibit 82 the FDA was
19 granting Arthrex the permission to sell Arthrex
20 FiberWIRE; is that right?

21 A Correct.

22 Q And Exhibit 82 is dated May 14th, 2001; is that
23 right?

24 A Correct.

25 Q Did Arthrex sell FiberWIRE prior to May 14th.

BROOKSTEIN DECLARATION EXHIBIT 19



April 26, 2001.

Food and Drug Administration
Center for Devices and Radiological Health
Attn: Mr. David Krause
Department of Health and Human Services, Public Health Service
Division of General & Restorative Devices
9200 Corporate Boulevard
Rockville, Maryland 20850

**RE: Amendment to Original Pre-Market Submission 510(k) #K010673,
Arthrex FiberWIRE™**

On February 28, 2001 the Arthrex FiberWIRE™ was submitted to the FDA. Per a conversation between Ann Waterhouse (Arthrex) and David Krause (FDA) on April 16, 2001, the following is being sent as amendment information. Specifically included are percentages for content, labeling amendments, and research data for accessory equipment.

Arthrex has submitted in duplicate the requested information and respectfully asks that these be accorded the same confidentiality as the original submission, K010673. We request that the Food and Drug Administration keep confidential all information outside of the 510(k) summary and indications for use.

Should the following information be in any way deficient, please let us know. We will be happy to provide you with any missing details or information. Should you have further questions concerning the amendment we have submitted, please contact either Vernon Brown or Ann Waterhouse at (941) 643-5553. Thank you.

Sincerely,

Ann Waterhouse
Regulatory Affairs Specialist

CONFIDENTIAL - NON-PATENT
PROSECUTION COUNSEL
ONLY

ARM 002103

1. In many places in the document the suture is described as being composed of polyester and ultrahigh molecular weight polyethylene (UHMWPE). Please describe the polyester and indicate what percentage of the suture weight is UHMWPE. Also, what percentage of the weight of the suture is the dye and what percentage of the weight of the suture is the coating?

Description of the polyester:

The Polyester used by Pearsalls to produce Arthrex FiberWIRE™ is created from high tenacity filaments of Polyethylene Terephthalate. Specifically, this is type 712 polyester that is manufactured by KoSa GmbH & Co KG.

% of the suture weight which is Ultra high molecular weight polyethylene:

Broken into percentage, the Polyester is 38.09% of the suture input and the Polyethylene (UHMW) is 61.91%.

Suture weight which is dye:

In using an accepted dye, D&C Blue No. 6, neither Pearsalls nor Arthrex measured the percentage of weight which the suture gained by the dye process. Pearsalls certifies that the process with which they dye the polyester conforms with the 2.0 Cupric Sulfate listed in the USP Matching Solutions table on page 1585 of USP 24. Also, Pearsalls certifies to all of it's customers that the D&C Blue No. 6 is a FDA approved dye.

Suture weight which is coating:

Pearsalls makes the following statement about the NuSil coating applied to both our Arthrex FiberWIRE™ and other competitor product:

"The coating is a silicone rubber known as NuSil 2045. It is identical to the coating we apply for Davis & Geck silk suture, and is also used to coat "Ticron" polyester suture. As noted (above) we cannot measure the amount of coating so the product is accepted by our customers on the basis of an agreed detailed coating process which includes mixing the NuSil to a certain viscosity and the speeds, temperatures, and other parameters for the coating process. For each coating batch all these details are recorded in the batch documentation available to you and other customers."

CONFIDENTIAL - NON-PATENT
PROSECUTION COUNSEL
ONLY

000001

ARM 002104

2. You describe your suture as having a silicone elastomer coating. Please identify a legally marketed suture predicate that is coated with silicone.

The predicate devices with Nu-Sil silicone derivative or equivalent, used to coat silk and polyester suture are listed below. The 510 (k) summaries or statements for these predicate products are contained in the following pages:

K930586	Dermalon, Surgilon, Ophthalon, & Ophthalmic Suture
K930590	Silk Sutures
K961925	Polyester Non-Absorbable Surgical Suture
K990088	Synthofil Non-Absorbable PET Surgical Suture
K001172	Polyester Non-Absorbable Surgical Suture
K003590	Grams Polyester Non-Absorbable Suture

CONFIDENTIAL - NON-PATENT
PROSECUTION COUNSEL
ONLY


ARM 002105

000002

BROOKSTEIN DECLARATION EXHIBIT 20



BROOKSTEIN DECLARATION EXHIBIT 21



15KV

48x

208µm

0000